
Bair Hugger*

Temperature Management Unit

Model 750

Operator's Manual

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
Total Temperature Management System

Bair Hugger

Temperature Management Unit
Model 750

Operator's Manual

English	5
Deutsch	21
Français	39
Italiano	57
Español	75
Português	93
Nederlands	111
Dansk	129
Svenska	147
Norsk	165
Sumoi	183
Ελληνικά	201
日本語	219
Polski	237
Magyar	255
Česky	273
Lietuvių Kalba	291

Operator's Manual

5

English

Table of Contents

Introduction.....	7
Description of the Total Temperature Management System.....	7
Indications	7
Contraindications	7
Warnings	7
Cautions	8
Proper Use and Maintenance	9
Read Before Servicing Equipment.....	9
Overview and Operation.....	10
Unit Power-On Reset.....	10
Overview of Control Panel.....	11
Mounting the Temperature Management Unit to an IV Pole	13
Instructions for Use	14
Viewing the Temperature Mode Timer.....	15
What to Do in Case of an Over-Temperature Condition	15
What to Do in Case of a Fault Condition	15
General Maintenance.....	16
Cleaning the Cabinet and Hose	16
Technical Support and Customer Service.....	17
When You Call for Technical Support	17
In-Warranty Repair and Exchange	17
Specifications	18
Definition of Symbols	20

Introduction

Description of the Total Temperature Management* System

The Bair Hugger brand Total Temperature Management system consists of a Model 750 forced-air temperature management unit (also available with rolling stand and sheet clip) and disposable components, including Bair Hugger forced-air blankets, Bair Paws* patient warming gowns and the 241* blood/fluid warming set. You can use the Model 750 temperature management unit in all clinical settings including the operating room to provide patient temperature management.

This manual includes operating instructions and unit specifications for the Model 750 temperature management unit. For information about using Bair Hugger blankets, Bair Paws gowns or the 241 blood/fluid warming set with the Model 750 temperature management unit, refer to the “Instructions for Use” included with each of these disposable components.

Indications

The Bair Hugger Temperature Management system is intended to prevent and treat hypothermia. In addition, the temperature management system can be used to provide patient thermal comfort when conditions exist that may cause patients to become too warm or too cold. The temperature management system can be used with adult and pediatric patients.

Contraindications

Do not apply heat to lower extremities during aortic cross-clamping. Thermal injury may occur if heat is applied to ischemic limbs.

Warnings

1. Do not leave patients with poor perfusion unmonitored during prolonged warming therapy sessions. Thermal injury may result.
2. The Bair Hugger temperature management unit has been designed to operate safely ONLY with Arizant Healthcare* disposable components. Use with other products may cause thermal injury. (To the full extent permitted by law, the manufacturer and/or importer declines all responsibility for thermal injury resulting from the unit being used in conjunction with non-Arizant Healthcare products.)
3.  Do not warm patients with the temperature management unit hose alone. Thermal injury may result. Always connect the hose to a Bair Hugger blanket or Bair Paws gown before providing therapy.
4. Do not place the non-perforated side of the blanket on the patient. Thermal injury may result. Always place the perforated side (that is, with small holes) towards the patient.
5. Do not continue temperature management therapy if the over-temp indicator light illuminates and the alarm sounds. Thermal injury may result. Unplug the unit, and contact a qualified service technician.

6. Do not continue 241° blood/fluid warming therapy if the over-temp indicator light illuminates and the alarm sounds. Immediately stop fluid flow, and discard the blood/fluid warming set. Unplug the temperature management unit, and contact a qualified service technician.
7. Do not use a forced-air warming device over transdermal medication. Increased drug delivery and patient injury or death may occur.
8. Do not allow the patient to lie on the warming unit hose or allow the hose to directly contact the patient's skin during patient warming; thermal injury may result.
9. Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
10. Resuable blankets made from woven fabric, or blankets without discrete, visible holes, can cause the safety system of this unit to fail, which may result in serious thermal injury. This warming unit has been designed to operate safety ONLY with Bair Hugger® blankets and Bair Paws gowns.

Cautions

1. Except for specific blanket models, Bair Hugger blankets are not sterile and are all intended for single patient use ONLY. Placing a sheet between the Bair Hugger blanket and the patient does not prevent contamination of this product.
2. Monitor the patient's temperature and cutaneous response every 10-20 minutes or according to institutional protocol and monitor the patient's vital signs regularly. Adjust air temperature or discontinue therapy when the therapeutic goal is reached or if vital sign instability occurs. Notify physician of vital sign instability immediately.
3. Do not leave pediatric patients unattended during therapy.
4. Do not initiate temperature management therapy unless the temperature management unit is safely placed on a hard surface or securely mounted. Otherwise, injury may result.
5. To prevent tipping, clamp the temperature management unit to an IV pole at a height that provides stability. We recommend clamping the unit no higher than 44" (112cm) on an IV pole with a minimum 14" (35cm) radius wheelbase. Failure to do so may result in IV pole tipping, catheter site trauma, and patient injury.
6. Electrical shock hazard. Do not disassemble the temperature management unit unless you are a qualified service technician. There are electrically live parts within the unit when it is connected to a power source, even when the unit is in Standby mode.

Notices

1. The Bair Hugger temperature management unit meets medical electronic interference requirements. If radio frequency interference with other equipment should occur, connect the unit to a different power source.
2. Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare professional.
3. To reliably ground this Bair Hugger warming unit, only connect to receptacles marked "Hospital Only" or "Hospital Grade".

Proper Use and Maintenance

Arizant Healthcare Inc. assumes no responsibility for the reliability, performance, or safety of the temperature management unit if the following events occur:

- Modifications or repairs are performed by unqualified personnel.
- The unit is used in a manner other than that described in the Operator's or Service Manual.
- The unit is installed in an environment that does not meet the appropriate electrical and grounding requirements.

Read Before Servicing Equipment

All repair, calibration, and servicing of the temperature management unit require the skill of a qualified, medical equipment service technician who is familiar with good practice for medical device repair. If service does not require the manufacturer's attention, the Model 750 Service Manual provides the technical information needed to service the unit, or Arizant Healthcare Inc. will provide that information on request. Perform all repairs and maintenance in accordance with the instructions in the Service Manual.

Overview and Operation


Figure A. Model 750 temperature management unit – front view

Figure B. Model 750 temperature management unit – rear view

Unit Power-On-Reset

The Model 750 temperature management unit proceeds through the following power-on-reset sequence after you connect the unit to a grounded power source:

- Performs all self-test functions.
- Illuminates indicator lights and all pixels in the alphanumeric display momentarily.
- Displays the text “BH 750” and the software revision level in the alphanumeric display.
- Sounds the alarm (three low-level clicks).
- Enters the *Standby* mode.

If the temperature management unit loses power for less than 1 second, the unit's software will return the unit to the operating mode you selected prior to the power loss. However, if the temperature management unit loses power for longer than 1 second, the unit's software will reset when you restore power. The unit will then enter the *Standby* mode.

Overview of Control Panel


Figure C. Model 750 temperature management unit control panel

ALPHANUMERIC DISPLAY

The alphanumeric display shows the temperature at the blanket end of the temperature management unit hose in degrees Celsius.

TEMPERATURE MODES

- Press the Low (32°C), Med (38°C), or High (43°C) button to select the desired temperature.
- Press the Ambient button to supply room temperature air.
- When you select a temperature mode, the following events occur:
 - Corresponding indicator light illuminates.
 - Blower operates.
 - Temperature at the blanket end of the hose appears in the alphanumeric display.
 - Heater activates except in Ambient mode.
 - The timer resets itself every time you select a different temperature mode.

STANDBY

When the unit is in Low, Medium, High, or Ambient mode, press the Standby button to place the unit in *Standby* mode.

When you select the Standby mode, the following events occur:

- Standby indicator light illuminates.
- Blower and heater are turned OFF.
- Alphanumeric display deactivates.
- Alarm and fault detection functions remain active.
- Operating timer pauses.

TEMP IN RANGE INDICATOR LIGHT

The Temp in Range indicator light illuminates when the temperature at the blanket end of the hose is $\pm 1.5^{\circ}\text{C}$ of the selected setting; this indicator light does not illuminate in the Ambient mode.

FAULT INDICATOR LIGHT

When a system fault occurs, the amber Fault indicator light flashes and an alarm sounds.

Please refer to *What to Do in Case of a Fault Condition* on page 15 for additional information.

OVER-TEMP INDICATOR LIGHT

If the unit senses an over-temperature condition, the red *Over-temp* indicator light flashes and an alarm sounds.

Please refer to *What to Do in Case of an Over-Temperature Condition* on page 15 for additional information.

Mounting the Temperature Management Unit to an IV Pole

CAUTION

To prevent tipping, clamp the Model 750 temperature management unit to an IV pole at a height that provides stability. We recommend clamping the unit no higher than 44 in. (112 cm) from the floor on an IV pole with a minimum 14 in. (35.6 cm) radius wheelbase. Failure to do so may result in IV pole tipping, catheter site trauma, and patient injury.

METHOD

To mount the Model 750 temperature management unit to an IV pole, simply turn the clamp handle clockwise to tighten the clamp on the IV pole and counterclockwise to release it. See Figure D.


Figure D.

Instructions for Use

The instructions below describe how to operate the Model 750 temperature management unit. For information about using Bair Hugger blankets or the 241 blood/fluid warming set with the Model 750 temperature management unit, refer to the Instructions for Use included with each of these disposable components.

1. If the Model 750 temperature management unit is not mounted on an IV pole or Bair Hugger rolling stand, place the unit on a flat, hard surface, such as a table, before beginning temperature management therapy. Do not place the temperature management unit on a soft or uneven surface, such as a bed, or the air intake may become blocked and cause the unit to overheat.

2. Insert the end of the temperature management unit's hose into the hose port of a Bair Hugger blanket. Use a twisting motion to ensure a snug fit.

3. Connect the temperature management unit to a properly grounded power source. The unit will be operating in Standby mode and the Standby indicator light will illuminate.
4. Press the appropriate button to select the desired temperature. When the unit reaches the selected temperature, the Temp in Range indicator light will illuminate; this indicator light does not illuminate in the Ambient mode.
5. Monitor the patient's temperature and cutaneous response every 10-20 minutes or according to institutional protocol and monitor the patient's vital signs regularly. Adjust air temperature or discontinue therapy when the therapeutic goal is reached or if vital sign instability occurs. Notify physician of vital sign instability immediately.
6. Press the Standby button when temperature management therapy is complete and discard disposable components.
7. Unplug the temperature management unit to disconnect it from the power source.

Viewing the Temperature Mode Timer

The temperature mode timer records the amount of time the temperature management unit has been operating in the selected temperature mode. The timer resets itself every time you select a different temperature mode.

To view the temperature mode timer, press the button for the current temperature mode, and hold it for 2 seconds. The alphanumeric display will show the temperature mode time, then go back to showing the temperature mode setting.

What to Do in Case of an Over-Temperature Condition

If an over-temperature condition occurs, the red *Over-temp* indicator light flashes and an alarm sounds. The unit automatically turns the heater, blower, and operating mode indicator lights OFF. The control panel will not respond to commands, but pressing any button will silence the alarm.

IF AN OVER-TEMPERATURE CONDITION OCCURS:

1. Discontinue all temperature management therapy. If you are using the 241 blood/fluid warming system, immediately stop fluid flow and discard the blood/fluid warming set.
2. Unplug the temperature management unit.
3. Contact a qualified service technician.

What to Do in Case of a Fault Condition

The Model 750 temperature management unit's software recognizes several non-hazardous conditions and reports those conditions as faults. When a system fault occurs, the software stores the fault code in memory, the amber Fault indicator light flashes, and an alarm sounds. The control panel will not respond to commands, but pressing any button will silence the alarm.

IF A FAULT CONDITION OCCURS:

1. Unplug the temperature management unit and wait for 5 minutes.
2. Reconnect the temperature management unit to a grounded power source. The unit will perform the normal power-on-reset sequence and then enter the *Standby* mode.
3. Reselect the temperature setting.
4. If the unit does not return to normal operation, contact a qualified service technician.

General Maintenance

Cleaning the Cabinet and Hose

WARNING

Do not immerse the cabinet or hose while cleaning. Moisture will damage the components, and thermal injury may result.

PRECAUTIONS

- Do not use a dripping wet cloth to clean the cabinet. Moisture may seep into the electrical contacts and damage the components.
- Do not use alcohol or other solvents to clean the cabinet. Solvents may damage the labels and other plastic parts.

METHOD

1. Disconnect the temperature management unit from the power source before cleaning.
2. Wipe the cabinet and the outside of the hose with a damp, soft cloth and a mild detergent or antimicrobial spray.
3. Dry with a separate soft cloth.

Technical Support and Customer Service

USA, WORLDWIDE

TEL:

800-733-7775

952-947-1200

FAX:

800-775-0002

952-947-1400

WITHIN EUROPE

TEL:

+49-4154-9934-0

0800-100-1324 (Toll-free in Germany)

0800-877-077 (Toll-free in Switzerland)

FAX:

+49-4154-9934-20

0800-100-1324 (Toll-free in Germany)

0800-877-088 (Toll-free in Switzerland)

When You Call for Technical Support

We will need to know the serial number of your Bair Hugger temperature management unit when you call us. The serial number label is located on the back or side of the Bair Hugger temperature management unit.

In-Warranty Repair and Exchange

USA

Call Arizant Healthcare customer service if your Model 750 temperature management unit requires factory service. A customer service representative will give you a Return Authorization (RA) number. Please use this RA number on all correspondence concerning your temperature management unit. Your customer service representative will also send a shipping carton to you at no charge, if needed. Call your local supplier or sales representative to inquire about borrowing a temperature management unit while we service your unit. For more detailed instructions on returning units for service, please refer to the Model 750 Service Manual.

WORLDWIDE

Contact your local distributor concerning in-warranty repair and exchange.

Specifications

Physical Characteristics

DIMENSIONS	12.5" high x 13.5" deep x 10" wide 31.5 cm high x 34 cm deep x 25 cm wide
WEIGHT	16.3 lb (7.4 kg)
RELATIVE NOISE LEVEL	55 dBA
HOSE	Flexible, compatible with the Bair Hugger brand 241 blood/fluid warming system.
FILTRATION SYSTEM	0.2 µm air filter
RECOMMENDED FILTER CHANGE	Every 6 months or 500 hours of use.
MOUNTING	Can be clamped to an IV pole, placed on a hard surface, or mounted to the rolling stand accessory.

Temperature Characteristics

RECOMMENDED OPERATING ENVIRONMENT	15°C-25°C
TEMPERATURE CONTROL	Electronically controlled.
HEAT GENERATED	1644 BTU/hr (average), 482 W (average)
OPERATING TEMPERATURES	Average temperatures at the end of the hose: HIGH: 43° ± 1.5°C 109.4° ± 2.7°F MED: 38° ± 1.5°C 100.4° ± 2.7°F LOW: 32° ± 1.5°C 89.6° ± 2.7°F

Safety System

THERMOSTAT	Independent electronic circuit; thermal cut-off shuts the heater OFF at preset high temperature of 53°C ± 3°C at the end of the hose, back-up over-temperature detection at hose inlet.
ALARM SYSTEM	Over-temperature (53° ± 3°C): red Over-temp indicator light flashes, alarm sounds, heater and blower shut down, operating indicator lights turn OFF, control panel becomes unresponsive. Fault: amber Fault indicator light flashes, alarm sounds.
OVERCURRENT PROTECTION	Dual input fused lines.

Electrical Characteristics

HEATING ELEMENT	1400 W Resistive
LEAKAGE CURRENT	Meets UL 60601-1 and IEC 60601-1 requirements.
BLOWER MOTOR	Operating speed: approximately 4,000 rpm Airflow: up to 48 cfm or 23 L/s
POWER CONSUMPTION	Peak: 1550 W Average: 800 W
POWER CORD	15 ft., SJT, 3 cond., 13 A 15 ft., SJT, 3 cond., 15 A 4.6 m, HAR, 3 cond., 10 A
DEVICE RATINGS	110-120 VAC, 50/60 Hz, 11.7 A, or 220-240 VAC, 50/60 Hz, 7.2 A, or 100 VAC, 50/60 Hz, 15 A
FUSES	12 A (110 - 120 VAC) 8 A (220 - 240 VAC) 15 A (100 VAC)
CERTIFICATIONS	IEC 60601-1; EN 60601-1-2; UL 60601-1; CAN/CSA-C22.2, No. 601.1, EN 55011; ASTM F2196-02
CLASSIFICATION	Classified under IEC 60601-1 Guidelines (and other national versions of the Guidelines) as Class I, Type BF, Ordinary equipment, Continuous operation. Not suitable for use in the presence of flammable anesthetic mixtures with air or with oxygen or nitrous oxide. Classified by Underwriters Laboratories Inc. with respect to electric shock, fire, and mechanical hazards only, in accordance with UL 60601-1, ASTM F2196-02 and in accordance with Canadian/CSA C22.2, No. 601.1. Classified under the Medical Device as a Class IIb device.
DIAGNOSTICS	A qualified service technician can perform over-temperature detection system testing, temperature output testing, operating temperature calibration, and fault code troubleshooting.

Definition of Symbols

The following symbols may appear on the product's labeling or exterior packaging.

	ON/STANDBY
	ON (used on isolation switch)
	OFF (used on isolation switch)
	ON/OFF push button switch
	Temperature control
	Equipotentiality plug (ground)
	Fuse
	Caution Attention, consult appropriate documents
	Warning
	Nonexplosion proof
	Dangerous voltage
	Type BF equipment (patient applied)
	Voltage, alternating current (AC)
	Special Refuse, discard separately
	Protective earth ground
	No free hosing

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Authorized Representative in the European Community:
Arizant UK Ltd., Calder Island Way, Wakefield, WF2 7AW, United Kingdom
TEL (44) 1924 200550 • FAX (44) 1924 200518

Bair Hugger temperature management products may be covered by one or more of the following patents: US 7,101,389; 6,355,915; 6,309,409; 6,309,408; 6,290,716; 6,254,337; 6,241,755; 6,228,107; 6,203,567; 6,176,870; 6,168,612; 6,146,412; 6,129,936; 6,126,681; 6,126,393; 6,036,722; 5,997,572; 5,964,792; 5,928,274; 5,824,025; 5,800,489; 5,773,275; 5,733,318; 5,697,963; 5,674,269; 5,658,325; 5,620,482; 5,545,194; RE38462; 5,350,417; 5,336,250; 5,324,320; 5,300,102; 5,300,101; 5,184,612; 5,044,361; Japan 2,561,326; Canada 1,325,484; Australia 756,900; EPO 0311336; Other patents pending.

*Bair Hugger, Arizant, Arizant Healthcare, bright ideas that work, and the Bair Hugger and Arizant logos are trademarks of Arizant Healthcare Inc., registered or pending in the U.S. Patent & Trademark Office and in other countries.

Arizant
Healthcare
bright ideas that work

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*Les noms Bair Hugger, 241 et Total Temperature Management, et les logos Bair Hugger et Arizant sont des marques d'Arizant Healthcare Inc., déposées et en cours d'homologation à l'Office américain des brevets et des marques (U.S. Patent and Trademark Office) et dans d'autres pays.
*Bair Hugger, Arizant, bright ideas that work, Bair Hugger and Arizant logos, 241, and Total Temperature Management are trademarks of Arizant Healthcare Inc., registered or pending in the U.S. Patent & Trademark Office and in other countries.